

Model 20/50 Concentrate Disinfectant

The 20/50 supports a quaternary disinfectant. The disinfectant is mixed with four gallons of water in the main sink. To meet the required dilution ratio, 2 oz of disinfectant will be dispensed into the main sink. This is done automatically by the system.

In addition to the dilution ratio, a typical contact time of 5 to 10 minutes is required. Contact time is also done by the system.

The use of our No-Foam is occasionally needed based on the hardness or softness of the local water supply.

To ensure the disinfectant is being applied, a flow sensor monitors flow of the disinfectant liquid. If there is no flow of the fluid, the system will alert the operator and stop. A replacement of the empty container needs to be done by the operator.

A Safety Data Sheet is available on request.